


## Anti-Human Cd134/Ox40 Monoclonal Antibody (Pe) For Flow Cytometry - P43489

Item Number: CM0012560



### Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD134/OX40, validated for flow cytometry. This costimulatory receptor is key for T-cell immunity research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human CD134/OX40 mAb                                                                                     |
| Remarks/Alias    | OX40; ACT35; CD134; IMD16; TXGP1L                                                                                       |
| Species          | Human                                                                                                                   |
| GeneID           | 7293                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 29-216 of human CD134/OX40 (NP_003318.1). |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | FC                                                                                                                      |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | P43489                                                                                                                  |
| Protein Weight   | 29kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |