

Anti-Human Cd138 Pe Monoclonal Antibody For Flow Cytometry - P18827

Item Number: CM0011867



Introduction

Rabbit monoclonal antibody against human CD138/Syndecan-1, PE-conjugated, suitable for flow cytometry (FC). Designed for research use in cancer and immunology studies involving exosome biology and cell surface proteoglycan analysis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD138 mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan
Species	Human
Gene ID (Human)	6382
Gene ID	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-251 of human CD138/Syndecan-1 (NP_002988.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag