

Anti-Human Cd138/Syndecan-1 Monoclonal Antibody (Conjugate) For Flow Cytometry - P18827

Item Number: CM0009925



Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human CD138/Syndecan-1 (SWISS: P18827). Validated for flow cytometry (FC), suitable for studying cell surface proteoglycans in exosome and membrane biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD138/Syndecan-1 mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan
Species	Human
Gene ID (Human)	6382
Gene ID	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-251 of human CD138/Syndecan-1 (NP_002988.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag