

## Anti-Human Cd138/Syndecan-1 Rabbit Monoclonal Antibody For Flow Cytometry - P18827

Item Number: CM0012725



### Introduction

-conjugated rabbit anti-human CD138/Syndecan-1 monoclonal antibody for direct flow cytometry detection. Targets a cell-surface proteoglycan involved in exosome biogenesis and cell-matrix interactions. Validated on human samples.

[Learn More](#)

| Parameter        | Value                                                                                                                        |
|------------------|------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD138/Syndecan-1 mAb                                                                                       |
| Remarks/Alias    | SDC, CD138, SYND1, syndecan                                                                                                  |
| Species          | Human                                                                                                                        |
| GeneID           | 6382 (Human)                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 23-251 of human CD138/Syndecan-1 (NP_002988.4) |
| Source           | Rabbit                                                                                                                       |
| Category         | Monoclonal Antibody                                                                                                          |
| Application      | FC (Flow Cytometry)                                                                                                          |
| Cross Reactivity | Human                                                                                                                        |
| SWISS            | P18827                                                                                                                       |
| Protein Weight   | 32 kDa                                                                                                                       |
| Shipping         | Ice bag                                                                                                                      |