


Anti-Human Cd140B/Pdgfr β Monoclonal Antibody For Flow Cytometry - P09619

Item Number: CM0005939



Introduction

Rabbit monoclonal antibody against human CD140b/PDGFR β , validated for flow cytometry. A research tool for studying PDGF receptor signaling, vascular development, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD140b/PDGFR β mAb
Remarks/Alias	PDGFRB; CD140B; IBGC4; IMF1; JTK12; KOG5; PDGFR; PDGFR-1; PDGFR1; PENTT; platelet-derived growth factor receptor beta
Species	Human
GeneID	5159
Immunogen	Recombinant protein - Recombinant fusion protein containing a sequence corresponding to amino acids 33-530 of human CD140b/PDGFR β (NP_002600.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	P09619
Protein Weight	37 kDa / 123 kDa
Shipping	Ice bag