

Anti-Human Cd146 Monoclonal Antibody () For Flow Cytometry - P43121

Item Number: CM0009829



Introduction

Rabbit monoclonal CD146 antibody conjugated to for flow cytometry detection of human CD146/MUC18 in endothelial and melanoma research. Valid for IVD raw material supply.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD146 mAb
Remarks/Alias	CD146; MUC18; HEMCAM; METCAM; MelCAM
Species	Human
Gene ID	4162
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag