

Anti-Human Cd150/Slam Monoclonal Antibody () For Fc - Q13291

Item Number: CM0012259



Introduction

-conjugated rabbit monoclonal antibody targeting human CD150/SLAM (SLAMF1) for flow cytometry. Suitable for studying SLAM-mediated immune signaling, T/B cell activation, and measles virus receptor interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD150/SLAM mAb
Remarks/Alias	SLAM; CD150; CDw150
Species	Human
GeneID	6504
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-237 of human CD150 (NP_003028.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13291
Protein Weight	37kDa
Shipping	Ice bag