

Anti-Human Cd152/Ctla-4 Monoclonal Antibody For Fc (Intra) - P16410

Item Number: CM0005796



Introduction

High-quality rabbit anti-human CD152/CTLA-4 monoclonal antibody for intracellular flow cytometry (FC intra). Targets the inhibitory T-cell receptor CTLA-4, validated for human samples. Ideal for IVD and immune-oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD152/CTLA-4 mAb
Remarks/Alias	CD; GSE; GRD4; ALP55; CD152; CTLA-4; IDDM12; CELIAC3
Species	Human
GeneID (Human)	1493
GeneID	1493
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 37-162 of human CD152/CTLA-4 (NP_005205.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P16410
Protein Weight	6kDa/8kDa/19kDa/24kDa
Shipping	Ice bag