

Anti-Human Cd153 Monoclonal Antibody For Flow Cytometry - P32971

Item Number: CM0010742

Introduction

-conjugated rabbit monoclonal antibody recognizing human CD153 (CD30 ligand). Validated for flow cytometry; suitable for studying CD30/CD153-mediated T-cell proliferation and immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD153 mAb
Remarks/Alias	CD153; CD30L; CD30LG; TNLG3A
Species	Human
Gene ID (Human)	944
Gene ID	944
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 63-234 of human CD153 (NP_001235.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32971
Protein Weight	26kDa
Shipping	Ice bag