

Anti-Human Cd158 (Kir2DL1/S1/S2/S3/S5) Monoclonal Antibody For Flow Cytometry - Q14954 / Q14953 / Q14952 / P43631 / P43626

Item Number: CM0010052

Introduction

Fluorophore-conjugated rabbit monoclonal antibody targeting human CD158 (KIR2DL1/S1/S2/S3/S5). Validated for flow cytometry, it enables NK cell receptor profiling and immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD158 (KIR2DL1/S1/S2/S3/S5) mAb
Remarks/Alias	NKAT9; CD158G□ CD158A□ NKAT1; CD158H; CD158J; NKAT5
Species	Human
Gene ID (Human)	3810/3808/3806/3802/100132285
Gene ID	3810/3808/3806/3802/100132285
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human CD158 (NP_055327.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
SWISS	Q14954/Q14953/Q14952/P43631/P43626
Protein Weight	34kDa
Shipping	Ice bag