

## Anti-Human Cd158A/Kir2DL1 Monoclonal Antibody ( ) For Flow Cytometry - P43626

Item Number: CM0010660



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD158a/KIR2DL1 (p58.1, NKAT1). Validated for flow cytometry (FC). Used in NK cell inhibition, immune checkpoint, and receptor signaling studies.

[Learn More](#)

| Parameter        | Value                                                                                                                      |
|------------------|----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD158a/KIR2DL1 mAb                                                                                       |
| Remarks/Alias    | NKAT; NKAT1; p58.1; CD158A; KIR221; NKAT-1; KIR-K64; KIR2DL3                                                               |
| Species          | Human                                                                                                                      |
| Gene ID          | 3802                                                                                                                       |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human CD158a/KIR2DL1 (NP_055033.2) |
| Source           | Rabbit                                                                                                                     |
| Category         | Monoclonal Antibodies                                                                                                      |
| Application      | FC                                                                                                                         |
| Cross Reactivity | Human                                                                                                                      |
| SWISS            | P43626                                                                                                                     |
| Protein Weight   | 38kDa/41kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                                    |