

Anti-Human Cd15/Ssea-1 Monoclonal Antibody For Flow Cytometry - P22083

Item Number: CM0008097

Introduction

Mouse anti-Human CD15/SSEA-1 monoclonal antibody for flow cytometry. Recognizes Lewis X antigen (CD15) on human cells; protein weight ~59 kDa. Suitable for leukocyte immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD15/SSEA-1 mAb
Remarks/Alias	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV
Species	Human
Gene ID (Human)	2526
Gene ID	2526
Immunogen	Synthetic peptide A synthesized peptide derived from human CD15/SSEA-1.
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22083
Protein Weight	59 kDa
Shipping	Ice bag