


Anti-Human Cd166/Alcam Monoclonal Antibody For Flow Cytometry - Q13740

Item Number: CM0009061



Introduction

-conjugated rabbit monoclonal antibody targeting human CD166/ALCAM (SWISS: Q13740), validated for flow cytometry. Ideal for immunological studies involving T-cell activation, cell adhesion, and hematopoiesis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD166/ALCAM mAb
Remarks/Alias	MEMD; CD166
Species	Human
GeneID (Human)	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-526 of human CD166/ALCAM(NP_001618.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13740
Protein Weight	65kDa
Shipping	Ice bag