


Anti-Human Cd166/Alcam Monoclonal Antibody For Flow Cytometry - Q13740

Item Number: CM0011733



Introduction

Rabbit anti-Human CD166/ALCAM monoclonal antibody validated for Flow Cytometry (FC). Reacts with Human. Recombinant immunogen. Recognizes ALCAM/CD166, a marker for T-cell activation and stem cell research. Ideal for immune cell adhesion and oncology studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD166/ALCAM mAb
Remarks/Alias	MEMD; CD166
Species	Human
GeneID (Human)	214
GeneID	214
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-526 of human CD166/ALCAM (NP_001618.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13740
Protein Weight	15kDa/33kDa/63kDa/65kDa
Shipping	Ice bag