


Anti-Human Cd170/Siglec-5 Monoclonal Antibody For Fc - O15389

Item Number: CM0011323



Introduction

Rabbit monoclonal antibody targeting human CD170/Siglec-5, conjugated with ABflo 594, for flow cytometry. Binds specifically to CD170 on myeloid cells, enabling immunophenotyping and sialic acid binding studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD170/Siglec-5 mAb
Remarks/Alias	CD170; OBBP2; CD33L2; OB-BP2; SIGLEC-5
Species	Human
GeneID (Human)	8778
GeneID	8778
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 17-434 of human CD170/Siglec-5 (NP_003821.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O15389
Protein Weight	61kDa
Shipping	Ice bag