


Anti-Human Cd172B (Sirp-Beta-1) Monoclonal Antibody For Flow Cytometry - 000241

Item Number: CM0010049



Introduction

-conjugated rabbit anti-human CD172b monoclonal antibody, specific for human CD172b, validated for flow cytometry (FC). Suitable for studies of myeloid cell activation and negative regulation of receptor tyrosine kinase-coupled signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD172b mAb
Remarks/Alias	CD172b; SIRP-BETA-1
Species	Human
Gene ID (Human)	10326
Gene ID	10326
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-371 of human CD172b (NP_006056.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O00241
Protein Weight	43kDa
Shipping	Ice bag