

Anti-Human Cd177 Monoclonal Antibody For Flow Cytometry - Q8N6Q3

Item Number: CM0011199

Introduction

-conjugated rabbit monoclonal antibody targeting human CD177 (HNA-2a, PRV-1). Validated for flow cytometry (FC) with specificity for human. Ideal for neutrophil subset analysis in research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD177 mAb
Remarks/Alias	CD177; HNA-2a; HNA2A; NB1; NB1 GP; PRV-1; PRV1; NB1GP; CD177 antigen
Species	Human
GeneID (Human)	57126
GeneID	57126
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-407 of human CD177(NP_065139.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8N6Q3
Protein Weight	15kDa/27kDa/46kDa
Shipping	Ice bag