

Anti-Human Cd2 Monoclonal Antibody () For Flow Cytometry - P06729

Item Number: CM0010788



Introduction

-conjugated rabbit anti-human CD2 monoclonal antibody for flow cytometry. Detects T-cell surface antigen CD2 (LFA-2, T11) on T and NK cells. Swiss-Prot P06729.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD2 mAb
Remarks/Alias	T11; SRBC; LFA-2
Species	Human
Gene ID	914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag