

Anti-Human Cd2 Monoclonal Antibody () For Flow Cytometry - P06729

Item Number: CM0010716



Introduction

-conjugated rabbit monoclonal antibody against human CD2, optimized for flow cytometry. Ideal for studying T-cell adhesion, activation, and immunophenotyping. Validated for human reactivity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD2 mAb
Alias	T11, SRBC, LFA-2
Species	Human
Gene ID	914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-209 of human CD2 (NP_001758.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag