

Anti-Human Cd225/Ifitm1 Monoclonal Antibody For Flow Cytometry (Intracellular) - P13164

Item Number: CM0008203

Introduction

-conjugated rabbit anti-human CD225/IFITM1 monoclonal antibody validated for intracellular flow cytometry. This interferon-induced transmembrane protein restricts viral entry and is overexpressed in multiple cancers. Ideal for antiviral and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD225/IFITM1 mAb
Remarks/Alias	IFITM1; 9-27; CD225; DSPA2a; IFI17; LEU13; interferon-induced transmembrane protein 1
Species	human
GeneID	8519
Immunogen	Recombinant fusion protein containing a sequence within amino acids 1-80 of human CD225/IFITM1 (NP_003632.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Human
SWISS	P13164
Protein Weight	13kDa
Shipping	Ice bag