

Anti-Human Cd225/Ifitm1 Monoclonal Antibody () For If/Icc And Flow Cytometry - P13164

Item Number: CM0010811

Introduction

-conjugated rabbit monoclonal antibody against human CD225/IFITM1, validated for immunofluorescence and intracellular flow cytometry. Suitable for virology and oncology research; detects a 14kDa protein (P13164).

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD225/IFITM1 mAb
Remarks/Alias	9-27; CD225; IFI17; LEU13; DSPA2a
Species	Human
GeneID	8519
Immunogen	Recombinant fusion protein containing a sequence within amino acids 1-80 of human CD225/IFITM1 (NP_003632.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	P13164
Protein Weight	14kDa
Shipping	Ice bag