


## Anti-Human Cd226/Dnam-1 Monoclonal Antibody For Flow Cytometry - Q15762

Item Number: CM0004690



### Introduction

Rabbit monoclonal antibody targeting human CD226/DNAM-1 (DNAM-1), validated for flow cytometry. Immunogen spans extracellular domain (aa 19-247). Cross-reacts with human. Ideal for studying immune cell adhesion and cytotoxicity.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD226/DNAM-1 mAb                                                                                        |
| Remarks/Alias    | PTA1; DNAM1; DNAM-1; TLISA1                                                                                               |
| Species          | Human                                                                                                                     |
| GeneID (Human)   | 10666                                                                                                                     |
| GeneID           | 10666                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 19-247 of human CD226/DNAM-1 (NP_006557.2). |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | FC                                                                                                                        |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | Q15762                                                                                                                    |
| Protein Weight   | 39kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |