

Anti-Human Cd239/Bcam () Monoclonal Antibody For Flow Cytometry - P50895

Item Number: CM0012086



Introduction

Rabbit monoclonal antibody against human CD239/BCAM, conjugated to . Validated for flow cytometry. Target is a cell adhesion molecule with roles in migration and invasion.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD239/BCAM mAb
Remarks/Alias	BCAM
Species	Human
Gene ID (Human)	4059
Gene ID	4059
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-547 of human CD239/BCAM (NP_005572.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P50895
Protein Weight	67kDa
Shipping	Ice bag