

Anti-Human Cd24 Monoclonal Antibody () For If/Icc And Fc - P25063

Item Number: CM0009274

Introduction

Conjugated rabbit monoclonal antibody against human CD24, suitable for immunofluorescence and flow cytometry. Targets a GPI-anchored glycoprotein involved in B-cell activation, differentiation, and immune modulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD24 mAb
Remarks/Alias	CD24A; CD24
Species	Human
GeneID (human)	100133941
GeneID	100133941
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-80 of human CD24 (NP_037362.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P25063
Protein Weight	8kDa; 13kDa
Shipping	Ice bag