

## Anti-Human Cd244/2B4 Monoclonal Antibody For Flow Cytometry - Q9Bzw8

Item Number: CM0009759

### Introduction

Rabbit monoclonal anti-human CD244/2B4 antibody conjugated to . Designed for flow cytometry, this antibody targets the SLAM family receptor CD244 on NK and T cells. UniProt Q9BZW8.

[Learn More](#)

| Parameter        | Value                                                                                                          |
|------------------|----------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD244/2B4 mAb                                                                                |
| Remarks/Alias    | CD244; 2B4; NAIL; NKR2B4; Nmrk; SLAMF4; CD244 molecule                                                         |
| Species          | Human                                                                                                          |
| GenelD           | 51744                                                                                                          |
| Immunogen        | Recombinant protein corresponding to a sequence within amino acids 22-221 of human CD244/2B4 (NP_001160135.1). |
| Source           | Rabbit                                                                                                         |
| Category         | Monoclonal Antibodies                                                                                          |
| Application      | Flow Cytometry                                                                                                 |
| Cross Reactivity | Human                                                                                                          |
| SWISS            | Q9BZW8                                                                                                         |
| Protein Weight   | 30kDa/36kDa/41kDa                                                                                              |
| Shipping         | Ice bag                                                                                                        |