

Anti-Human Cd27 Monoclonal Antibody For Flow Cytometry - P26842

Item Number: CM0007791



Introduction

-conjugated rabbit anti-human CD27 monoclonal antibody for flow cytometry (FC) detection of human CD27, a T-cell activation antigen and immune checkpoint receptor. Targets extracellular domain (aa 21-192). 29kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD27 mAb
Remarks/Alias	T14; S152; Tp55; TNFRSF7; S152. LPFS2
Species	Human
GeneID (Human)	939
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-192 of human CD27 (NP_001233.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P26842
Protein Weight	29kDa
Shipping	Ice bag