


Anti-Human Cd271/Ngfr Monoclonal Antibody For Fc - P08138

Item Number: CM0008067



Introduction

Rabbit monoclonal antibody conjugated to , specifically recognizing human CD271/NGFR (p75NTR). Validated for flow cytometry, suitable for detecting neurotrophin receptor expression in neuronal and immune cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD271/NGFR mAb
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR
Species	Human
GeneID (Human)	4804
GeneID	4804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-250 of human CD271/NGFR (NP_002498.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08138
Protein Weight	35kDa/45kDa
Shipping	Ice bag