

Anti-Human Cd3 (Cd3E) Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0004108

Introduction

Mouse anti-Human CD3 monoclonal antibody for flow cytometry (FC). Targets the human CD3 epsilon chain (CD3E), a key subunit of the T-cell receptor-CD3 complex essential for adaptive immunity. Suitable for human samples. P07766.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GenelD (Human)	916
GenelD	916
Immunogen	Recombinant protein of human CD3
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag