

Anti-Human Cd3 Epsilon Monoclonal Antibody For Flow Cytometry - P07766

Item Number: CM0008985



Introduction

Rabbit anti-Human CD3 epsilon monoclonal antibody validated for flow cytometry. Targets human CD3 epsilon chain (23kDa), part of TCR-CD3 complex essential for T-cell activation and development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD3 mAb
Remarks/Alias	T3E; TCRE; IMD18; CD3epsilon
Species	Human
GeneID (Human)	916
GeneID	916
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-126 of human CD3 (NP_000724.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07766
Protein Weight	23kDa
Shipping	Ice bag