


Anti-Human Cd30 (Tnfrsf8) Monoclonal Antibody () For Flow Cytometry - P28908

Item Number: CM0009999



Introduction

Rabbit monoclonal anti-human CD30 (TNFRSF8) antibody conjugated to . Validated for flow cytometry to detect the CD30 receptor involved in lymphoblast growth and transformation. Ideal for immune oncology and lymphoma research.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD30 mAb
Remarks/Alias	TNFRSF8; CD30; DIS166E; Ki-1; TNF receptor superfamily member 8
Species	Human
GeneID	943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-379 of human CD30 (NP_001234.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28908
Protein Weight	14kDa/50kDa/63kDa
Shipping	Ice bag