


Anti-Human Cd300F/Irem-1 Monoclonal Antibody For Fc - Q8Tdq1

Item Number: CM0004341



Introduction

Rabbit monoclonal antibody against human CD300F/IREM-1, validated for flow cytometry. Recognizes the inhibitory receptor CMRF35-like molecule 1, involved in immune homeostasis and efferocytosis. Ideal for immunology research.

[Learn More](#)

Property	Value
Product Name	Rabbit anti-Human CD300F/IREM-1 mAb
Remarks/Alias	CD300LF; CD300f; CLM-1; CLM1; IREM-1; IREM1; IgSF13; LMIR3; NKIR; CMRF35-like molecule 1
Species	Human
Gene ID (Human)	146722
Gene ID	146722
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-156 of human CD300F/IREM-1 (NP_620587.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8TDQ1
Protein Weight	18kDa/21kDa/26kDa/32kDa
Shipping	Ice bag