

Anti-Human Cd303/Bdca-2 Monoclonal Antibody For Flow Cytometry - Q8Wtt0

Item Number: CM0012483



Introduction

Rabbit anti-Human CD303/BDCA-2 mAb for flow cytometry. Specifically detects C-type lectin domain family 4 member C (BDCA-2), a plasmacytoid dendritic cell marker. Optimal for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD303/BDCA-2 mAb
Remarks/Alias	DLEC; HECL; BDCA2; CD303; BDCA-2; CLECSF7; CLECSF11; PRO34150
Species	Human
GeneID	170482 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 45-213 of human CD303/BDCA-2 (NP_569708.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8WTT0
Protein Weight	21kDa/25kDa
Shipping	Ice bag