


Anti-Human Cd304/Neuropilin-1 Monoclonal Antibody For Fc - O14786

Item Number: CM0010685



Introduction

-conjugated rabbit monoclonal antibody targeting human CD304/Neuropilin-1, validated for flow cytometry. Functions in angiogenesis, axon guidance, and as SARS-CoV-2 host factor. Ideal for developmental biology and infection research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD304/Neuropilin-1 mAb
Remarks/Alias	NRP1; BDCA4; CD304; NP1; NRP; VEGF165R; neuropilin-1
Species	Human
GeneID (Human)	8829
GeneID	8829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-644 of human CD304/Neuropilin-1 (NP_003864.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O14786
Protein Weight	68kDa/71kDa/103kDa
Shipping	Ice bag