

Anti-Human Cd31/Pecam1 Monoclonal Antibody For Fc - P16284

Item Number: CM0004335



Introduction

-conjugated rabbit monoclonal antibody targeting human CD31/PECAM1, validated for flow cytometry. Ideal for studying leukocyte transendothelial migration, endothelial cell adhesion, and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD31/PECAM1 mAb
Remarks/Alias	CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM
Species	Human
GeneID (Human)	5175
GeneID	5175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-601 of human CD31 (NP_000433.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16284
Protein Weight	83kDa
Shipping	Ice bag