


Anti-Human Cd32 Monoclonal Antibody For Flow Cytometry - P12318

Item Number: CM0012611



Introduction

-conjugated rabbit monoclonal antibody targeting human CD32 (FCGR2A). Suitable for flow cytometry detection of low affinity IgG Fc receptor on monocytes, neutrophils, and platelets. Recognizes extracellular domain.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD32 mAb
Remarks/Alias	CD32; FCG2; FcGR; CD32A; CDw32; FCGR2; IGFR2; FCGR2A1; FcgammaRIIa
Species	Human
Gene ID (Human)	2212
Gene ID	2212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-218 of human CD32 (NP_001129691.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12318
Protein Weight	35kDa
Shipping	Ice bag