

Anti-Human Cd337/Nkp30 Monoclonal Antibody For If/Icc, Fc - O14931

Item Number: CM0008516



Introduction

Rabbit monoclonal antibody against human CD337/NKp30 (NCR3), validated for IF/ICC and FC (intra) applications. Recognizes the natural cytotoxicity triggering receptor expressed on NK cells, involved in anti-tumor immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD337/NKp30 mAb
Remarks/Alias	NCR3; 1C7; CD337; LY117; MALS; NKp30; natural cytotoxicity triggering receptor 3
Species	Human
Gene ID	259197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of human CD337/NKp30 (NP_667341.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	O14931
Protein Weight	16kDa/17kDa/18kDa/19kDa/20kDa/21kDa
Shipping	Ice bag