

Anti-Human Cd34 Monoclonal Antibody () For If/Icc And Fc - P28906

Item Number: CM0010253

Introduction

Rabbit monoclonal antibody against human CD34, directly conjugated to . Suitable for immunofluorescence and flow cytometry. Ideal for research on hematopoietic progenitor cells and vascular endothelium markers.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Human
GeneID (Human)	947
GeneID	947
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	35kDa/41kDa
Shipping	Ice bag