

Anti-Human Cd34 Monoclonal Antibody (Pe/Cyanine7) For Flow Cytometry - P28906

Item Number: CM0004908

Introduction

PE/Cyanine7-conjugated rabbit monoclonal antibody against human CD34, validated for flow cytometry. Targets the hematopoietic progenitor cell antigen, a 41 kDa membrane glycoprotein involved in cell adhesion and stem cell homing.

[Learn More](#)

Parameter	Value
Product Name	PE/Cyanine7 Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3
Species	Human
Gene ID	947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	41kDa
Shipping	Ice bag