


## Anti-Human Cd352/Slamf6 Monoclonal Antibody For Flow Cytometry - Q96Du3

Item Number: CM0004756



### Introduction

Rabbit monoclonal -conjugated antibody targeting human CD352/SLAMF6 (SLAM family member 6). Validated for flow cytometry, enabling detection of NTB-A on NK, T, and B cells for immune response studies.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD352/SLAMF6 mAb                                                                                          |
| Remarks/Alias    | SLAMF6; CD352; KALI; KALib; Ly108; NTB-A; NTBA; SF2000; SLAM family member 6                                                |
| Species          | Human                                                                                                                       |
| GeneID           | 114836 (Human)                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 24-215 of human CD352/SLAMF6 (NP_001171643.1) |
| Source           | Rabbit                                                                                                                      |
| Category         | Monoclonal Antibodies                                                                                                       |
| Application      | FC                                                                                                                          |
| Cross Reactivity | Human                                                                                                                       |
| SWISS            | Q96DU3                                                                                                                      |
| Protein Weight   | 24kDa/37kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                                     |