


Anti-Human Cd355/Crtam Monoclonal Antibody For Flow Cytometry - 095727

Item Number: CM0010651



Introduction

Rabbit anti-Human CD355/CRTAM monoclonal antibody for flow cytometry. CD355/CRTAM mediates heterophilic cell adhesion regulating T-cell and NK cell activation, cytotoxicity, and tissue retention.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD355/CRTAM mAb
Remarks/Alias	CD355
Species	Human
GeneID (Human)	56253
GeneID	56253
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-287 of human CD355/CRTAM (NP_062550.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	095727
Protein Weight	22kDa/44kDa
Shipping	Ice bag