


Anti-Human Cd38 Monoclonal Antibody For Fc - P28907

Item Number: CM0009500



Introduction

Rabbit anti-Human CD38 monoclonal antibody for flow cytometry (FC) applications, validated for human samples. Ideal for studying ADP-ribosyl cyclase/cyclic ADP-ribose hydrolase 1 (CD38) in research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag