

Anti-Human Cd39/Entpd1 Monoclonal Antibody For Flow Cytometry - P49961

Item Number: CM0005088

Introduction

Rabbit monoclonal antibody conjugated to for flow cytometric detection of human CD39/ENTPD1, a membrane ectonucleotidase regulating nucleotide signaling, expressed on activated lymphocytes and endothelial tissues.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD39/ENTPD1 mAb
Remarks/Alias	CD39; SPG64; ATPDase; NTPDase-1
Species	Human
GeneID (Human)	953
GeneID	953
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 255-442 of human CD39/ENTPD1 (NP_001767.3) .
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P49961
Protein Weight	58kDa
Shipping	Ice bag