

Anti-Human Cd45Ra Polyclonal Antibody For If-P - P08575

Item Number: CM0021684



Introduction

Rabbit polyclonal antibody targeting human CD45RA (PTPRC), validated for immunofluorescence-paraffin (IF-P). Suitable for immune cell typing and T-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Human CD45RA Rabbit pAb
Remarks/Alias	PTPRC; B220; CD45; CD45R; GP180; L-CA; LCA; LY5; T200; receptor-type tyrosine-protein phosphatase C; Human CD45RA
Species	Human
GeneID (human)	5788
GeneID	5788
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 10-110 of human CD45RA(P08575-8).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P
Cross Reactivity	Human
SWISS	P08575
Protein Weight	130kDa/147kDa
Shipping	Ice bag