


Anti-Human Cd5 Monoclonal Antibody For Flow Cytometry - P06127

Item Number: CM0011052



Introduction

Rabbit anti-human CD5 monoclonal antibody validated for flow cytometry. Targets T-cell surface glycoprotein CD5 (T1/Leu-1), a negative regulator of TCR signaling with roles in immune tolerance. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD5 mAb
Remarks/Alias	T1; LEU1
Species	Human
Gene ID (human)	921
Gene ID	921
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag