


Anti-Human Cd55 (Daf) Rabbit Monoclonal Antibody () For Flow Cytometry - P08174

Item Number: CM0005275



Introduction

Rabbit monoclonal anti-human CD55 antibody conjugated with . Validated for flow cytometry on human samples. Targets complement decay-accelerating factor (DAF), a key regulator of complement activation. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD55 mAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE
Species	Human
Gene ID (Human)	1604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-353 of human CD55 (NP_000565.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08174
Protein Weight	39kDa/41kDa/48kDa/56kDa/59kDa
Shipping	Ice bag