

Anti-Human Cd56 Monoclonal Antibody For Flow Cytometry - P13591

Item Number: CM0009001



Introduction

Mouse monoclonal antibody against human CD56 (NCAM1), validated for flow cytometry. Relevant for neuroscience and viral infection research, including neuron adhesion and rabies/Zika receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Mouse anti-Human CD56 mAb
Remarks/Alias	CD56; NCAM; MSK39
Species	Human
GeneID (Human)	4684
GeneID	4684
Immunogen	Recombinant protein of human CD56
Source	Mouse
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13591
Protein Weight	40kDa/73kDa/80kDa/83kDa/93kDa/95kDa
Shipping	Ice bag