


## Anti-Human Cd59 Monoclonal Antibody For Flow Cytometry - P13987

Item Number: CM0010695



### Introduction

-conjugated rabbit anti-human CD59 monoclonal antibody for flow cytometry. Targets CD59 (protectin, MAC-inhibitory protein), a key complement regulator, for immune research applications.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD59 mAb                                                                                             |
| Remarks/Alias    | 1F5; EJ16; EJ30; EL32; G344; MIN1; MIN2; MIN3; M1RL; HRF20; MACIF; MEM43; MIC11; MSK21; 16.3A5; HRF-20; MAC-IP; p18-20 |
| Species          | Human                                                                                                                  |
| Gene ID (Human)  | 966                                                                                                                    |
| Gene ID          | 966                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 26-101 of human CD59 (NP_000602.1)       |
| Source           | Rabbit                                                                                                                 |
| Category         | Monoclonal Antibodies                                                                                                  |
| Application      | FC                                                                                                                     |
| Cross Reactivity | Human                                                                                                                  |
| SWISS            | P13987                                                                                                                 |
| Protein Weight   | 14kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |