

Anti-Human Cd62L/L-Selectin Monoclonal Antibody For Fc - P14151

Item Number: CM0010844



Introduction

-conjugated rabbit monoclonal antibody against human CD62L/L-selectin. Validated for flow cytometry. Designed for studying leukocyte adhesion, lymphocyte homing, and inflammation. - P14151.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD62L/L-Selectin mAb
Remarks/Alias	SELL; CD62L; LAM1; LECAM1; LEU8; LNHR; LSEL; LYAM1; PLNHR; TQ1; L-selectin
Species	Human
GeneID (Human)	6402
GeneID	6402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14151
Protein Weight	42kDa/43kDa
Shipping	Ice bag