

## Anti-Human Cd68 Monoclonal Antibody (Percp) For Flow Cytometry (Intracellular) - P34810

Item Number: CM0012804

### Introduction

PerCP-conjugated rabbit monoclonal antibody against human CD68 (macrosialin, gp110). Suitable for intracellular flow cytometry detection of a 37 kDa lysosomal/endosomal membrane protein in macrophages and tumor cells.

[Learn More](#)

| Parameter        | Value                                                                                                                                 |
|------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PerCP Rabbit anti-Human CD68 mAb                                                                                                      |
| Remarks/Alias    | GP110; LAMP4; SCARD1                                                                                                                  |
| Species          | Human                                                                                                                                 |
| GeneID (Human)   | 968                                                                                                                                   |
| GeneID           | 968                                                                                                                                   |
| Immunogen        | Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 22-319 of human CD68(NP_001242.2)) |
| Source           | Rabbit                                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                                 |
| Application      | FC (intra)                                                                                                                            |
| Cross Reactivity | Human                                                                                                                                 |
| SWISS            | P34810                                                                                                                                |
| Protein Weight   | 37kDa                                                                                                                                 |
| Shipping         | Ice bag                                                                                                                               |