

Anti-Human Cd7 Monoclonal Antibody (Biotin) For Flow Cytometry - P09564

Item Number: CM0011518

Introduction

Biotinylated rabbit monoclonal antibody targeting human CD7, validated for flow cytometry (FC). Recognizes a 25 kDa transmembrane glycoprotein on T-cells and NK cells. Suitable for immunophenotyping and T-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Biotin Rabbit anti-Human CD7 mAb
Remarks/Alias	GP40; TP41; Tp40; LEU-9
Species	Human
GeneID (Human)	924
GeneID	924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag